**DPP® Syphilis Screen & Confirm Assay**  
**FOR EXPORT ONLY**  
FOR *IN VITRO* DIAGNOSTIC USE  
FOR PROFESSIONAL USE ONLY

A Qualitative Test for the Simultaneous Detection of Treponemal and Non-Treponemal antibodies in Human Fingertick Whole Blood, Venous Whole Blood, Serum, or Plasma Specimens.

Read this Product Insert completely before using the product. Follow the instructions carefully when performing the test as not doing so may result in inaccurate Test Results.

STORAGE: Store at 2 to 30°C (36 to 86°F)

#### INTENDED USE

The Chembio DPP® Syphilis Screen & Confirm Assay is a single-use immunochromatographic rapid screening test for the simultaneous detection of antibodies against non-Treponemal and *Treponema pallidum* antigens in fingertick whole blood, venous whole blood, serum, and plasma. The Chembio DPP Syphilis Screen & Confirm Assay is intended for use as a point-of-care (POC) test to aid in the screening and confirmation of infection with Syphilis.

#### SUMMARY AND EXPLANATION

Syphilis is a sexually transmitted disease (STD) caused by the spirochete *Treponema pallidum*. It is a chronic bacterial infection that remains a public health concern worldwide, especially in resource poor settings. Syphilis can be transmitted from infected women to their unborn children during pregnancy. Worldwide 12 million individuals are diagnosed with syphilis each year, 90% of them in developing countries. Diagnosed individuals are also at risk of being coinfecting with and transmitting the HIV virus[1].

Early and appropriate diagnosis and treatment prevents the transmission and development of severe complications. A rapid serologic test for specific antibodies to non-Treponemal and *T. pallidum* antigens is important in the early diagnosis and treatment monitoring of syphilis patients. In turn, this monitoring allows for the formulation of a more successful public health strategy. Various serologic tests are currently available such as Venereal Disease Research Laboratory (VDRL), rapid plasma reagin (RPR), fluorescent Treponemal antibody absorption (FTA-ABS) test, *T. pallidum* hemagglutination (TPHA) test, immunoenzymatic assay (EIA), *Treponema pallidum* particle agglutination (TPPA) test and Western blot (WB) test [2– 5]. The Chembio DPP® Syphilis Screen & Confirm Assay is a unique non-Treponemal and Treponemal rapid point-of-care test, which is simple and easy to use. The DPP® Syphilis Screen & Confirm Assay is a qualitative immunoassay for the detection of antibodies to the non-Treponemal and *T. pallidum* antigens in serum, plasma and whole blood.

#### PRINCIPLE OF THE TEST

The Chembio DPP® Syphilis Screen & Confirm Assay employs a unique combination of protein A and anti-human IgM antibody, which are conjugated to colloidal gold particles. It also utilizes a recombinant antigen of *T. pallidum* and synthetic antigens for non-Treponema, separately bound to the membrane solid phase. The sample is applied to the Well 1 along with the Running Buffer. Five minutes after adding the sample + Running Buffer to Well 1, additional Running Buffer is added to BUFFER Well 2. The Running Buffer allows gold conjugates to migrate and bind to the antibody/antigen complex at the test lines. If the antibodies to syphilis (Treponemal and non-Treponemal) are present in the sample then two pink/purple lines (1 and 2) are produced. In the absence of anti-Treponemal and anti-non-Treponemal antibodies, there are no pink/purple lines in the TEST (1 and 2) area. Unbound conjugated gold particles continue to migrate along the membrane and produce a pink/purple line in the CONTROL (C) area. This procedural control serves to demonstrate that the reagents have been properly applied and have migrated through the device.

#### MATERIALS PROVIDED

Each kit contains the items to perform 20 tests:

- 20 DPP® Syphilis Screen & Confirm Individually Pouched Test Devices
- 20 Disposable Microsafe® Tubes
- 1 DPP® Syphilis Running Buffer Bottle (6mL) – Red Cap
- 1 Product Insert

#### ACCESSORIES AVAILABLE

- Chembio DPP® Syphilis Screen & Confirm Assay Controls  
Each package contains:
  - 1 DPP® Syphilis Reactive Control
  - 1 DPP® Syphilis Nonreactive Control
  - 1 Product Insert

#### MATERIALS REQUIRED BUT NOT PROVIDED

- Clock, watch or other timing device
- Pipettor capable of delivering 5µL of sample for serum or plasma specimens
- Pipettor capable of delivering 10µL of sample for venous whole blood samples
- Sterile gauze (for fingertick samples only)
- Disposable gloves
- Antiseptic wipes
- Biohazard disposal containers
- Sterile Safety Lancet (for fingertick whole blood samples)
- Collection devices for samples other than fingertick

## WARNINGS

### For *IN VITRO* diagnostic use

1. Read the Product Insert completely before using this assay. Follow the instructions carefully as not doing so may result in inaccurate test results.
2. Use of this test kit with sample types other than those specifically approved for use with this device may result in inaccurate test results.
3. This test should be performed at 18 to 30°C (64 to 86°F). If stored refrigerated, ensure that the pouch is brought to operating temperature before performing testing.
4. If the test kit is stored at temperatures outside the storage temperature 2 to 30°C (36 to 86°F), or used outside the operating temperature 18 to 30°C (64 to 86°F), use the Chembio DPP Syphilis Screen & Confirm Assay Controls to ensure proper performance of the test.

## PRECAUTIONS

### Safety Precautions

1. Handle the samples and materials contacting samples as if capable of transmitting infection.
2. Do not eat, drink or smoke in the area where samples and kit reagents are handled. Avoid any contact between hands, eyes or mouth during sample collection and testing.
3. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when handling patient samples.
4. Dispose of all samples and materials used in the test procedure in a biohazard waste container. Lancets should be placed in a puncture-resistant container prior to disposal. The recommended method of disposal of biohazard waste is autoclaving for a minimum of 1 hour at 121°C. Disposable materials may be incinerated. Liquid wastes may be mixed with appropriate chemical disinfectants. A freshly prepared solution of 10% bleach (0.5% solution of sodium hypochlorite) is recommended. Allow 60 minutes for effective decontamination.  
**NOTE: Do not autoclave solutions that contain bleach.**
5. Use 10% bleach or other appropriate disinfectants to wipe all spills. The bleach solution should be made fresh each day.
6. For additional information refer to: Centers for Disease Control (CDC): Updated U.S. Public Health Service Guidelines for the Management of Occupational Exposures to HBV, HCV and HIV and Recommendations for Postexposure Prophylaxis and Centers for Disease Control (CDC): Updated U.S. Public Health Service Guidelines for the Management of Occupational Exposures to HIV and Recommendations for Postexposure Prophylaxis [6,7].

### Handling Precautions

1. If Desiccant Packet is missing, DO NOT USE. Discard test device and use a new test device.
2. Do not open the sealed foil pouch until just prior to use.
3. Do not use any test device if the pouch has been perforated.
4. Each test device is for single use only.
5. Do not use the test beyond the expiration date printed on the pouch. Always check expiration date prior to testing.
6. Do not mix reagents from different lot numbers of kits.
7. Adequate lighting is required to read the test results.

## STORAGE AND STABILITY

The DPP<sup>®</sup> Syphilis Screen & Confirm Assay devices should be stored in unopened pouches at 2 to 30°C (36 to 86°F). Do not freeze. Do not open pouch until you are ready to perform a test. When stored as indicated, test devices are stable until the expiration date marked on the pouch. Running Buffer should be stored at 2 to 30°C (36 to 86°F) in its original bottle. Do not use beyond the indicated expiration date.

## SAMPLE COLLECTION

The Chembio DPP<sup>®</sup> Syphilis Screen & Confirm Assay can be performed on fingerstick whole blood, venous whole blood, serum, and plasma specimens.

### Fingerstick Whole Blood

Following laboratory procedure, clean the finger of the person being tested with an antiseptic wipe. Allow the finger to dry thoroughly or wipe dry with a sterile gauze pad. Using a sterile lancet, puncture the skin just off the center of the finger and wipe away the first drop of blood with sterile gauze. Avoid squeezing the fingertip to accelerate bleeding as this may dilute the blood with excess tissue fluid. Collect the sample from the second drop with a disposable 10µL Microsafe<sup>®</sup> Tube (provided). Follow test procedure instructions.

### Venous Whole Blood

Draw blood following laboratory procedure for obtaining venous blood. Depending on use, collect sample in a tube containing citrate, heparin or EDTA. Be sure the tube of blood is well mixed before sampling. Follow test procedure instructions.

### Serum or Plasma

Draw blood following laboratory procedure for obtaining serum or plasma samples. Collect serum samples in tubes that do not contain any anticoagulant (serum). Collect plasma samples in tubes containing citrate, heparin, or EDTA. Collect sample in a clean container following standard laboratory procedures. Be sure that the tube of serum or plasma is well mixed before sampling. Follow test procedure instructions.

Venous whole blood, serum and plasma specimens may be tested immediately after collection. If specimens are not tested immediately, refrigerate them at 2 to 8°C (36 to 46°F) following collection. These specimens should be tested within 3 days of collection, otherwise serum or plasma samples should be frozen at -20°C (-4°F) or colder. Allow refrigerated samples to reach room temperature and mix gently before testing. DO NOT FREEZE WHOLE BLOOD!

### Specimen Shipping

If specimens are to be shipped, they should be packed in compliance with regulations covering the transportation of etiologic agents. Venous whole blood, serum and plasma specimens should be shipped refrigerated with cold packs or wet ice.


**TEST PROCEDURE**

**Kit component preparation**

All components for the Chembio DPP® Syphilis Screen & Confirm Assay are ready to use as supplied. Follow directions as indicated. If the sample and / or kit components have been refrigerated, remove them from the refrigerator and allow them to come to a temperature of 18 to 30° C (64 to 86°F) prior to testing.


**TEST PROCEDURE**

1. Remove the Chembio DPP® Syphilis Screen & Confirm test device from its pouch and place it on a flat surface (it is not necessary to remove the desiccant from the pouch). **NOTE:** If desiccant packet is missing, DO NOT USE. Discard test device and use a new test device. Label the test device with patient ID or identification number (see Figure 1 below). Note that the DPP® test device has 2 colored lines in the Test Window; one is blue and the other is green. If either one of the 2 colored lines is absent, DO NOT USE. Discard the test device and use a new test device.

**2. Fingerstick Samples**

For fingerstick whole blood, prick the finger and wipe away the first drop. Collect the second drop with a Microsafe® Tube holding it in a horizontal position as shown. Touch the tip of the tube to the blood sample. Capillary action will draw the sample to the black fill line and stop.

**Never squeeze the bulb of the tube while sampling.**

**Venous Blood, Serum or Plasma Samples**

If testing a sample other than fingerstick blood, use an accurate laboratory pipette to obtain the required amount of sample.

1. Serum or plasma - 5µL
2. Whole blood - 10µL


3. Transfer the required amount of serum, plasma or whole blood into the center of the round SAMPLE + BUFFER Well 1 of the device.

For fingerstick blood, transfer the blood from the Microsafe Tube into the center of the round SAMPLE + BUFFER Well 1 of the device as shown in Figure 2 below.


**Figure 2**

To push the sample out of the Microsafe Tube, line up the tip of the tube with the SAMPLE + BUFFER Well 1 and squeeze the bulb.

**ONLY IF THE SAMPLE DOES NOT COME OUT OF THE TUBE,** hold the tube vertically and slide a finger over the vent hole near the black mark. Then line up the tip with the SAMPLE well and squeeze the bulb.


Well 1. Add sample here.

Well 1. Add sample here.

4. Invert the Running Buffer bottle (red cap), and hold it vertically (not at an angle) over SAMPLE + BUFFER Well 1. Add 2 drops of Running Buffer (~50 $\mu$ L) into the center of the round SAMPLE + BUFFER Well 1 (See Figure 3 below).  
**NOTE:** When removing the cap from the Buffer bottle, take care to remove only the red portion of the cap. Do not unscrew the white part from the bottle.


**Figure 3**

5. **Wait 5 minutes.** The blue and green colored lines should have disappeared from the rectangular window. If not, discard the test device and repeat the procedure with a new DPP test device.

Invert the Running Buffer bottle (red cap), and hold it vertically (not at an angle) over BUFFER Well 2. Slowly add 5 drops (~135 $\mu$ L) of Running Buffer to BUFFER Well 2 (see Figure 4 below). Start timing device.

**Figure 4**

6. Read the test result at 15 minutes after the addition of Running Buffer into BUFFER Well 2 (Step 4). In some cases a test line may appear in less than 15 minutes. However, 15 minutes are needed to report a non-reactive result. **Do not read results after 15 minutes from start of test (Start of test = First addition of Running Buffer added to Well 2).**

**NOTE:** Discard the used test device and any other test materials into a biohazard waste container.

#### LIMITATIONS OF THE PROCEDURE

1. The Chembio DPP<sup>®</sup> Syphilis Screen & Confirm Assay must be used with human capillary (fingerstick) or venous whole blood, serum, and plasma only. Using other types of specimens or testing of venipuncture whole blood specimens collected using a tube containing an anticoagulant other than citrate, heparin or EDTA may not yield accurate results. For serum samples, collect blood without anticoagulant.
2. The Chembio DPP<sup>®</sup> Syphilis Screen & Confirm Assay must be used in accordance with the instructions in this product insert to obtain accurate results.
3. Test results must be read at 15 minutes after the addition of Running Buffer to the BUFFER Well 2 (Step 5)
4. Ensure finger is completely dry before performing fingerstick.
5. Read results in a well-lit area.
6. A reactive result using the Chembio DPP<sup>®</sup> Syphilis Screen & Confirm Assay suggests the presence of antibodies to Syphilis *Treponema pallidum* and/or non-Treponemal antigens. The Chembio DPP<sup>®</sup> Syphilis Screen & Confirm Assay is intended as an aid in the diagnosis and treatment of Syphilis.
7. For a reactive result, the intensity of the test line does not necessarily correlate with the titer of antibody in the specimen.
8. A non-reactive result does not preclude the possibility of exposure to Syphilis or infection with Syphilis. An antibody response to recent exposure may take several months to reach detectable levels.

#### QUALITY CONTROL

##### Built-in Control Feature

The control line serves as a built-in internal control and gives confirmation of the presence of reagents and proper test performance. A pink/purple line will appear in the CONTROL (C) area if the test has been performed correctly and the device is working properly (Please see: Interpretation of Test Results).

**External Quality Control**





Chembio DPP<sup>®</sup> Syphilis Screen & Confirm Reactive and Nonreactive Controls are available separately for use with the Chembio DPP<sup>®</sup> Syphilis Screen & Confirm test. The Controls are used to verify the operator's ability to properly perform the test and to interpret the results. The Reactive Control will produce a reactive test result and has been manufactured to produce faint lines in the TEST (1)(2) area. The Nonreactive Control will produce a nonreactive test result. Run the controls as described in the Test Procedure section for a serum / plasma sample and follow the directions in the Interpretation of Results section of this product insert. It is the responsibility of each facility using the Chembio DPP<sup>®</sup> Syphilis Screen & Confirm Assay to establish an adequate quality assurance program to ensure the performance of the device under specific locations and conditions of use.

Run the DPP<sup>®</sup> Syphilis Screen & Confirm Assay Controls under the following circumstances:

- Each new operator prior to performing tests on patient samples
- When opening a new test kit lot
- Whenever a new shipment of test kits is received
- If the temperature of the test storage area falls outside of 2 to 30°C (36 to 86°F)
- If the temperature of the testing area falls outside of 18 to 30°C (64 to 86°F)
- At periodic intervals as indicated by the user facility

If the Control reagents do not produce the expected results, contact Chembio Diagnostic Customer Service at (1-800-327-3635).


**INTERPRETATION OF TEST RESULTS**

<p><b>NONREACTIVE:</b></p>	<p>One pink/purple line in the CONTROL (C) area with no line in the TEST (1) or TEST (2) areas indicates a nonreactive result. A nonreactive result at 15 minutes from the addition of Running Buffer to Well 2 indicates that there are no detectable antibodies to both Treponemal and non-Treponemal antigens. A nonreactive result does not exclude the possibility of infection with syphilis.</p>	 <p>No Test Lines      Control Line</p> <p>1   2   C</p>
<p><b>REACTIVE TREPONEMAL (Trep) and NON-TREPONEMAL (N.Trep):</b></p>	<p>Three pink/purple lines, one in the TEST (1) area, one in the TEST (2) area, and one in the CONTROL (C) area indicate a reactive result for Treponemal and non-Treponemal antibodies. The lines in the TEST areas may look different from the line in the CONTROL (C) area and the two TEST lines may look different from each other. Intensities of the TEST (1), TEST (2) and CONTROL (C) lines may vary. Test results with visible lines in TEST (1), TEST (2) and CONTROL (C) areas, regardless of intensity, are considered REACTIVE for both Treponemal and non-Treponemal antibodies. The test result is interpreted as Preliminary POSITIVE for Treponemal and non-Treponemal antibodies.</p>	 <p>(1) Trep Line   (2) N.Trep Line   Control Line</p> <p>1   2   C</p>
<p><b>REACTIVE TREPONEMAL (Trep) and NONREACTIVE NON-TREPONEMAL (N.Trep):</b></p>	<p>Two pink/purple lines, one in the TEST (1) area and one in the CONTROL (C) area indicate a reactive result for Treponemal antibodies. The line in the TEST area may look different from the line in the CONTROL (C) area. Intensities of the TEST (1) and CONTROL (C) lines may vary. Test results with visible lines in TEST (1) and CONTROL (C) areas, regardless of intensity, are considered REACTIVE for Treponemal antibodies. The test result is interpreted as Preliminary POSITIVE for Treponemal antibodies.</p>	 <p>(1) Trep Line      Control Line</p> <p>1   2   C</p>
<p><b>NONREACTIVE TREPONEMAL (Trep) and REACTIVE NON-TREPONEMAL (N.Trep):</b></p>	<p>Two pink/purple lines, one in the TEST (2) area and one in the CONTROL (C) area indicate a reactive result for non-Treponemal antibodies. The line in the TEST area may look different from the line in the CONTROL (C) area. Intensities of the TEST (2) and CONTROL (C) lines may vary. Test results with visible lines in TEST (2) and CONTROL (C) areas, regardless of intensity, are considered REACTIVE for non-Treponemal antibodies. The test result is interpreted as Preliminary POSITIVE for non-Treponemal antibodies.</p>	 <p>(2) N.Trep Line   Control Line</p> <p>1   2   C</p>

**INTERPRETATION OF TEST RESULTS (cont.)**

**INVALID:**

A pink/purple line should always appear in the CONTROL (C) area whether or not one or two lines appear in the TEST areas. If there is no distinct pink/purple line visible in the CONTROL (C) area, then the test is INVALID. An INVALID test cannot be interpreted. It is recommended that the INVALID test be repeated with a new device.

**EXPECTED RESULTS**

This is a qualitative test for the detection of antibodies to Treponemal and non-Treponemal antigens in whole blood, serum, and plasma. As described in the PERFORMANCE CHARACTERISTICS section below, the sensitivity of the Chembio DPP<sup>®</sup> Syphilis Screen & Confirm Assay was found to be substantially equivalent to RPR and TPPA tests when tested on selected performance panels.

**PERFORMANCE CHARACTERISTICS**

**Specificity**

Specificity was determined by evaluating the performance of the Chembio DPP<sup>®</sup> Syphilis Screen & Confirm Assay when tested with 90 known syphilis negative serum samples (Table 1). The Chembio DPP<sup>®</sup> Syphilis Screen and Confirm Assay had specificity of 99% (89/90) for the Treponemal and 100% (90/90) for the non-Treponemal band.

**Table 1**  
Specificity of DPP<sup>®</sup> Syphilis Screen & Confirm tested with 90 Syphilis Negative Sera

	# Samples Tested	# Positive	% Specificity
Treponemal	90	1	99%
Non-Treponemal	90	0	100%

**Performance with Brazilian Specimens**

The performance of the DPP<sup>®</sup> Syphilis Screen & Confirm Assay was evaluated with 35 positive and negative syphilis specimens obtained from Brazil and compared to the performance of two Treponemal assays (TPPA and EIA) and one non-Treponemal assay (RPR). These results are shown in Table 2 below.

**Table 2**  
Performance with 35 Brazilian Positive and Negative Specimens

Sample I.D.	Chembio DPP <sup>®</sup> Syphilis Screen & Confirm Lot: DPPSYP122910 Exp: 28 Oct 12		Fujirebio Seradia <sup>®</sup> TP-PA Lot: 19K09710 Exp: 6/30/2011	Trinity Captia <sup>™</sup> Syphilis IgG (s/co) Lot: 056 Exp: 12/31/2011	Fisher SureVue <sup>®</sup> RPR Lot: 0J02R3 Exp: 6/30/2012	
	Trep	Non Trep			Result	Titer
AEQ - 2C - 1	NR	NR	NR	0.39	Neg	--
AEQ - 2C - 2	R	R	R	3.22	Pos	4
AEQ - 2C - 3	R	R	R	4.94	Pos	32
AEQ - 2C - 4	R	R	R	2.79	Pos	8
AEQ - 2C - 5	NR	NR	NR	0.33	Pos	2
AEQ - 2C - 6	R	R	R	2.20	Pos	4
AEQ - 3C - 2	R	R	R	4.06	Pos	16
AEQ - 3C - 3	NR	NR	NR	0.37	Pos	2
AEQ - 3C - 4	R	R	R	2.57	Pos	16
AEQ - 3C - 5	NR	NR	NR	0.55	Neg	--
AEQ - 7C - 1	R	R	R	3.59	Pos	8
AEQ - 7C - 2	NR	NR	NR	0.35	Neg	--
AEQ - 7C - 3	R	R	R	3.18	Pos	16
AEQ - 7C - 4	NR	NR	NR	0.51	Neg	--
AEQ - 7C - 5	R	R	R	4.43	Pos	4
AEQ - 7C - 6	R	R	R	3.00	Pos	4
AEQ - 8C - 1	R	R	R	3.63	Pos	8
AEQ - 8C - 2	R	R	R	3.61	Pos	4
AEQ - 8C - 3	R	R	R	4.29	Pos	16
AEQ - 8C - 5	NR	NR	NR	0.53	Neg	--
AEQ - 8C - 6	R	R	R	2.92	Pos	2
AEQ - 9C - 1	NR	NR	NR	0.24	Pos	1
AEQ - 9C - 2	R	R	R	2.71	Pos	4
AEQ - 9C - 3	NR	NR	NR	0.44	Neg	--
AEQ - 9C - 4	R	R	R	3.00	Pos	8
AEQ - 9C - 5	NR	NR	NR	0.35	Neg	--
AEQ - 9C - 6	R	R	R	3.51	Pos	4
AEQ - 10C - 1	NR	NR	NR	0.37	Neg	--
AEQ - 10C - 2	NR	NR	NR	0.41	Neg	--
AEQ - 13C - 1	R	R	R	3.02	Pos	4
AEQ - 13C - 2	NR	NR	NR	0.37	Pos	1
AEQ - 13C - 3	NR	NR	NR	0.31	Neg	--
AEQ - 13C - 4	NR	NR	NR	0.41	Pos	1
AEQ - 13C - 5	R	R	R	3.37	Pos	4
AEQ - 13C - 6	NR	NR	NR	0.31	Neg	--

The DPP<sup>®</sup> Syphilis Screen & Confirm Assay produced identical results to the TPPA and EIA assays for the Treponemal line. The DPP<sup>®</sup> Screen & Confirm Assay agreed with the Fisher RPR assay with the exception of 5 samples (3 R1 and 2 R2) that were negative with the DPP<sup>®</sup> Syphilis Screen & Confirm Assay.

**Sensitivity - Mixed Titer Panel (BBI)**

The Chembio DPP<sup>®</sup> Syphilis Screen & Confirm Assay was tested against 20 members of Boston Biomedica, Inc. (BBI) mixed titer panel (PSS202) as shown in Table 3. The Chembio DPP<sup>®</sup> Syphilis Screen and Confirm Assay is similar to currently licensed assays in detecting these panel members.

**Table 3**

Comparison of Chembio DPP Syphilis Screen & Confirm Assay results with other licensed assays on 20 members of Seracare/BBI Mixed Titer Panel

Panel Member	Chembio DPP <sup>®</sup> Syphilis Screen & Confirm Test		Bect-Dick RPR	Wampole RPR	Diesse ATA	Fujirebio ATA	Olympus ATA	Trinity ATA
	Treponemal	Non-Treponemal	Titer	Titer	s/co	Result	Result	s/co
PSS202-01	Reactive	Reactive	128	128	>5.5	Reactive	Reactive	3.3
PSS202-02	Reactive	Reactive	8	8	>5.5	Reactive	Reactive	3.7
PSS202-03	Reactive	Reactive	4	4	>5.5	Reactive	Reactive	2.2
PSS202-04	Non-reactive	Non-reactive	Negative	Negative	0.3	Non-reactive	Non-reactive	0.2
PSS202-05	Reactive	Non-reactive	1	Negative	>5.5	Reactive	Reactive	2.25
PSS202-06	Reactive	Reactive	64	128	>5.5	Reactive	Reactive	4.6
PSS202-07	Reactive	Reactive	8	8	>5.5	Reactive	Reactive	3.5
PSS202-08	Reactive	Reactive	1	1	>5.5	Reactive	Reactive	0.2
PSS202-09	Reactive	Reactive	32	32	4.5	Reactive	Reactive	2.1
PSS202-10	Reactive	Reactive	16	64	>5.5	Reactive	Reactive	3.8
PSS202-11	Reactive	Reactive	4	4	5.1	Reactive	Reactive	3.1
PSS202-12	Reactive	Reactive	32	64	>5.5	Reactive	Reactive	2
PSS202-13	Reactive	Reactive	1	2	>5.5	Reactive	Reactive	3.1
PSS202-14	Reactive	Reactive	2	2	>5.5	Reactive	Reactive	3.7
PSS202-15	Reactive	Reactive	32	64	>5.5	Reactive	Reactive	2.2
PSS202-16	Non-reactive	Non-reactive	Negative	Negative	0.3	Non-reactive	Non-reactive	0.2
PSS202-17	Reactive	Reactive	2	2	>5.5	Reactive	Reactive	2.3
PSS202-18	Reactive	Non-reactive	1	Negative	>5.5	Reactive	Reactive	2.3
PSS202-19	Reactive	Reactive	2	2	>5.5	Reactive	Reactive	3.5
PSS202-20	Reactive	Reactive	Negative	Negative	>5.5	Reactive	Reactive	2.1
Run Date	11-Jan-11		29-Jan-03	19-Feb-03	27-Jan-03	14-Feb-03	28-Feb-03	4-Apr-03
Kit Lot #	SYP122910		225580	601	031-A	VN20905	UR20505	K21414A
Exp. Date	28-Oct-12		30-Sep-03	FEB 04	30-Jun-03	3-Sep-03	MAR 03	22-Aug-03
Product #	65-9520-0		4975539	6002	91100	6384	PH3000	800970Row

**PRECISION**
**Intraassay**

To be determined

**Interassay**

To be determined

**Cross Reactivity and Interference**

To be determined

**REFERENCES**

1. World Health Organization. (2001) Global prevalence and incidence of selected curable sexually transmitted infections overview and estimates. Geneva: WHO.
2. Sandra, A. L., Bret, M. S., and Andrew, H. R. (1995) Laboratory diagnosis and interpretation of tests for syphilis. Clin. Microbiol. Rev. 8(1):1-21.
3. Pope, V., M. B. Fears, W. E. Morrill, A. Castro, and S. E. Kikkert. (2000) Comparison of the Serodia *Treponema pallidum* particle agglutination, captia Syphilis-G, and SpiroTeck Reagin II tests with standard test techniques for diagnosis of syphilis. J. Clin. Microbiol. 38: 2543 - 2545.
4. Castro, R., E. S. Prieto, I. Santo, J. Azevedo, and F. da L. exposto. (2003) Evaluation of an enzyme immunoassay technique for detection of antibodies against *treponema pallidum*. J. Clin. Microbiol.41: 250 - 253.
5. Sambri, V., A. Marangoni, C. Eyer, C. Reichhuber, E. Soutschek, M. Negosanti, A. D'Antuono, and R. Cevenini. (2001) Western immunoblotting with five *Treponema pallidum* recombinant antigens for serological diagnosis of syphilis. Clin. Diagn. Lab. Immunol. 8:534 - 539.
6. Centers for Disease Control (CDC). (2001) Updated U.S. Public Health Service Guidelines for the Management of Occupational Exposures to HBV, HCV, and HIV and Recommendations for Postexposure Prophylaxis. MMWR. 50(RR11): 1-42.
7. Centers for Disease Control (CDC). (2005) Updated U.S. Public Health Service Guidelines for the Management of Occupational Exposure to HIV and Recommendations for Postexposure Prophylaxis. MMWR. 54(RR09): 1-17.

<b>SYMBOL LEGEND</b>	
	CONSULT THE MANUAL BEFORE USE
	CAUTION, CONSULT THE ACCOMPANYING DOCUMENTS
	DO NOT RE-USE
	FOR USE WITHIN TEMPERATURE LIMITS
	IN VITRO DIAGNOSTIC MEDICAL DEVICE
	BATCH CODE
	PRODUCT CATALOG NUMBER
	MANUFACTURERS IDENTIFICATION
	DATE OF MANUFACTURE
	USE BYE DATE

** CHEMBIO DIAGNOSTIC SYSTEMS, INC.  
 3661 HORSEBLOCK ROAD  
 MEDFORD, NY 11763 USA**

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