


DPP® (Dual Path Platform) Technology

LICENSING OPPORTUNITIES:


Ideal for Multiplexing Antigen and/or Antibody Detection

Independent delivery of test sample and conjugate leads to greater efficiency of specific immunodetection

Unique, Rapid, Point-of-Care Testing Platform



Improved visual interpretation, accuracy and reliability with “difficult” matrices (whole blood, oral fluid, etc.)




vs. Lateral Flow

Hand-held or phone-based reader (to “The Cloud”) provides objective interpretation and enables development of quantitative assays.

DPP® SampleTainer® specimen collection bottle:

A safe, closed system for handling potentially infectious blood and oral fluid samples.



U.S. and international patents issued and pending including U.S. patent 7,189,522

Please see chembio.com/innovation/patents

Contact: JEsfandiari@chembio.com

Phone: 631.924.1135 x112



DIAGNOSTIC SYSTEMS, INC.

Rapid Tests for Earlier Treatment™