

# HIV 1/2 STAT-PAK<sup>®</sup> Assay

A trusted, versatile point-of-care assay for the detection of HIV-1 and HIV-2 antibodies in fingerstick or venous whole blood, serum or plasma.

## Simple

- CLIA waived<sup>1</sup>
- 2-step test procedure
- Ready-to-use reagents
- Room temperature storage: 46° to 86°F (8° to 30°C)

## Reliable

- FDA approved
- Sensitivity: 99.7%; Specificity: 99.9%<sup>2</sup>
- Built-in IgG procedural control

## Fast

- Rapid results: 15 minutes

**Minimal sample size: 5 µL**



3661 Horseblock Road, Medford, NY 11763 USA  
1-844-CHEMBIO (1-844-243-6246)  
CustomerService@chembio.com

 **CHEMBIO**  
DIAGNOSTIC SYSTEMS, INC.  
*Rapid Tests for Earlier Treatment™*


<sup>1</sup> CLIA waived for fingerstick and venous whole blood  
<sup>2</sup> See product insert for complete performance details

**Manufactured in the USA**

## 2-Step Procedure\*



### 1 ADD SAMPLE

Transfer loop sample to Sample Well (S)



### 2 ADD RUNNING BUFFER

Add three drops of Running Buffer. Read results in 15 minutes.



\* Before testing, please refer to the complete product insert provided in the test kit

## Ordering Information

Reimbursement Code: 86703-QW, G0433-QW

Description	Contents	Catalog Number
HIV 1/2 STAT-PAK® Assay FDA Approved, CLIA Waived/Non-Waived	20 Individually Pouched Test Devices 20 Disposable 5µL Sample Loops 1 HIV Running Buffer (3.5 mL) 20 Copies of Subject Information Notice 1 Quick Reference Instructions (for CLIA Waived) 1 Product Insert	60-9505-1
HIV Rapid Test Control Pack FDA Approved	One Pack Contains: - 1 HIV-1 Reactive Control 0.25 mL - 1 HIV-2 Reactive Control 0.25 mL - 1 HIV Nonreactive Control 0.25 mL Each kit control vial contains sufficient volume to run 50 tests.	60-9549-0